

PAP.1 – EXTRA (Modified, Haematoxylin Solution Progressive Stain)

Code MS1457

Specifications

Description	A ready-to-use staining used for nucleus
	staining in cytology laboratory for the
	medical diagnostic purpose
Nuclei	Dark blue
Type of Staining	Progressive Stain
Mercury (Hg)	Free

Suitable for both manual staining and slide stainer for the Haematoxylin and Eosin stain (H&E stain)

TN_22-08-2021

S-QC03-MS-M-0103 / RV.00 / Eff.date 25-08-21